


PROTAGONISTA hypo-allergenic

Orthopaedic hypo-allergenic mattress with removable cover



Available in **minimalist** version, with no binding tape - as pages 40-41

This mattress is suitable for those who **are allergic to dust mites**, feel quite hot in bed but do not feel particularly cold, have moderate sweating, are overweight (exceed their ideal weight by 5/10%), **always sleep on their back and/or stomach**, do not wish to use an adjustable base, do not suffer from neck and/or back pain, do not have peripheral blood circulation issues, are not subject to electromagnetic pollution, want a mattress with differing comfort levels on either side (bi-comfort), have a particular need to remove the mattress cover.



Different **Aquapure® Soya** element firmness on either side of the mattress results in two comfort levels.



The **3D.V.S. fabric** perimetral insert facilitates ventilation inside the mattress.



The **bi-comfort** feature allows you to choose between two comfort levels.



AVAILABLE FIRMNESS LEVELS	
STANDARD	bi-comfort
	Medium-firm side (MeRi) Firm side (Ri)

STRUCTURE



Mini-bonnell springs
Innerspring structure made of 300 (on twin size) patented spring steel, quenched, hourglass mini-coils. Combined with **Aquapure® Soya**, with differing firmness on either side, it provides the mattress with both medium-firm or firm comfort levels (bi-comfort).



Adapto system
Anti-deformation and supporting material.



Isotex wear-resistant protection
Needle-punched, heat-set and sturdy.



Aquapure® Soya
components, convoluted surface with seven comfort zones and differing firmness on either side.

UNDERLINING



Cotton blend
Underlining in elasticized mixed cotton jersey.

COVER



Anti-dust mite fibres allergofiber®
Hypo-allergenic padding, treated with a permanent active anti-dust mites agent.



Cotton blend fabric Silverness® silver yarn cover in double elastic fabric with permanent anti-dust mite, anti-microbial anti-stress and anti-static effect.



3D.V.S. permeable band 3D Ventilation System
band on the mattress perimeter.




Removable cover
Four-sided sturdy zipper.



Dry-clean
The cover can be dry cleaned.



Washable at 60° C
The cover can be washed in water at 60° C.



Mattress height

28
303