


MICRO POCKETED SPRINGS AND MEMORY STD SOYA & ALOE SEVEN COMFORT ZONES



HEAD ZONE	SHOULDER ZONE	LOWER BACK ZONE	HIP ZONE		THIGH ZONE	CALF ZONE	FOOT ZONE
Intermediate support area, for the correct support to the head.	Softer area capable of effectively accommodate the shoulder and avoid blood circulation issues.	Important supporting and elastic area, capable of sustaining the lower back to encourage correct spine posture.	Controlled hip cushioning to correctly support the heaviest part of the body.		Elastic and supporting area, capable of sustaining the lower end of thighs.	Relatively soft area, suitable for supporting calves without compressing them, to prevent one of the main causes of tingling and swelling in the lower limbs.	Intermediate support level, for proper feet support.

The differing firmness levels on either side of the structure, in combination with **Memory std® Soya & Aloe** result in two comfort levels: soft and medium-soft (bi-comfort).



The exclusive **Memory std® Soya & Aloe** wave surface in seven variable patterns significantly increases comfort for the whole body. The empty spaces between coils, together with the shape of **Memory std® Soya & Aloe**, facilitate air flow inside the mattress.



The **Memory std® Soya & Aloe** element with exclusive wave surface and seven variable patterns provides the correct support of the entire body.



Memory std® Soya & Aloe is an innovative viscoelastic foam.



Controlled compression pockets

Each spring steel coil is barrel-shaped. Placed under pressure, i.e. compressed by about 20%, it is packed in durable, air permeable, and hypo-allergenic fabric.

A complex system of independent forces

Pocketed coil rows are then assembled with a process that creates perfect alignment, avoids contact, and ensures the individual movement of each single coil. A perimetral steel rod keeps the structure together.

Large number of coils for excellent performance

The higher the number of coils, the better performance and comfort. This simple rule is at the foundation of Somnium design for this technology. The user can enjoy the effect of over 610 micro-coils per square meter, the so-called 1000/2000 system.

Seven comfort zones

The use of coils with different gauges allows the mattress to respond dynamically to body shape and movement.



MICRO-POCKETED SPRINGS MATTRESSES WITH REMOVABLE COVER	
Ergonomic seasonal DOLCESONNO MICRO	Ergonomic hypo-allergenic LUNADORO MICRO