

MINI-BONNELL SPRINGS AND ELASTOCELL (massaging or standard)

Massaging version shown in the picture



HEAD ZONE	SHOULDER ZONE	LOWER BACK ZONE	HIP ZONE	THIGH ZONE	CALF ZONE	FOOT ZONE
<p>The best models in the Somnium orthopaedic collection use almost twice the number of coils used in standard spring mattresses: 300 for the single mattress ⁽⁵⁾ and 600 for the double mattress ⁽⁶⁾. The progressive adaptation of the springs, together with the increased support provided by their number, guarantees the benefits of firm support to the back and an anatomically correct posture.</p>				<p>The mattress supports the different weights of the body (back, hips, legs) in a balanced way, regardless of the total weight of the user. The density of the suspensions and their relatively small size is the result of an effort to increase density, which prevents wear of the insulating materials and paddings, extending the mattress life and significantly increasing comfort.</p>		
						<p>⁽⁵⁾ cm 85 x 190 ⁽⁶⁾ cm 165 x 190</p>

Comfort level: firm.



An effective relaxing and regenerating massage is made possible by the **Elastocell** convoluted surface with variable patterns. The empty spaces between coils, together with the convoluted surface (in the massaging version), facilitate air flow inside the mattress.

MINI-BONNELL SPRINGS AND ELASTOCELL (massaging or standard)

BONNELL SPRINGS MATTRESSES WITH NON-REMOVABLE COVER	
<p>Orthopaedic seasonal COMPLICE BRIO</p>	<p>Orthopaedic hypo-allergenic SMART ORTHOPEDIC</p>